



# KALOS

REPLACEMENT FLUID/BLOOD WARMER



FOR CRRT, HEMOPERFUSION, PLASMA TREATMENTS

Gruppo  
**MEDICA**  
Group

KALOS is easy to use and can be programmed by way of an operator panel equipped with:

- keyboard for setting the control set point to the desired temperature
- large backlit graphic display, allowing the simultaneous display of the outlet temperature, the set point value and information on the heater's functional status.
- acoustic alarms in the form of a buzzer and the display of explicit messages on the display.

**KALOS operates with a dedicated heating bag to be placed inside the device and connected to the extracorporeal circuits.**



Dimensions	260 x 145 x 80 mm
Weight	1.8 Kg
Power Supply	100 - 230 VAC selectable at the time of manufacture
Display	Liquid Crystal
Keyboard	4 Key Membrane
Heating element	200W.
Temperature of incoming liquid	[15 - 30] °C
Flow of incoming liquid	[10 - 200] ml/min
Temperature of the outlet fluid	Adjustable from 36 °C to 39 °C (Blood)
Accuracy	+/- 2 °C
Sensitivity	0.5 °C
Installation	Pole (150 ÷ 300 mm ext. diam)
Product Code:	M90117 - KALOS heating device 230V M90163 - KALOS heating Device 115V M90193 - KALOS Tubing set , ETO Sterile

M27369 rev 01 CDP050523B

Distribuito da:



MANUFACTURER

**MEDICA**

MEDICA S.p.A.  
Via degli Artigiani, 7  
41036, Medolla (Modena)  
Italia

Telefono +39 0535 51159  
Fax + 39 0535 52605  
info@medica.it  
www.medica.it



UNI EN ISO 9001: 2008  
UNI CEI EN ISO 13485: 2012

0476